

natus®



# Natus EP Supplies

Electrodes

## Reusable & Disposable Cup Electrodes

### Grass® Reusable EEG Stamped Cup Electrodes

- Choice of gold, silver or Ag/AgCl 10mm cups
- Pediatric 6mm Ag/AgCl cup available
- No-lift contour adjustable cup neck for low profile
- Molded plastic 1.5mm touchproof connector
- Six lead wire colors: two each: ●●●●●●

Part Number	Cup Type	Lead Length	Quantity
019-477500	10mm Gold	39" (1.0m)	12/pkg
019-477600	10mm Gold	59" (1.5m)	12/pkg
019-478600	10mm Gold	79" (2.0m)	12/pkg
019-477700	10mm Gold	98" (2.5m)	12/pkg
019-477800	10mm Silver	39" (1.0m)	12/pkg
019-477900	10mm Silver	59" (1.5m)	12/pkg
019-478700	10mm Silver	79" (2.0m)	12/pkg
019-478800	10mm Silver	98" (2.5m)	12/pkg
019-478400	10mm Ag/AgCl	39" (1.0m)	12/pkg
019-478500	10mm Ag/AgCl	59" (1.5m)	12/pkg
019-478900	10mm Ag/AgCl	79" (2.0m)	12/pkg
019-479000	10mm Ag/AgCl	98" (2.5m)	12/pkg
019-479100	6mm Ag/AgCl Pediatric	59" (1.5m)	12/pkg



### Reusable Snap Lead Pairs

- Designed for use with disposable snap electrodes
- Orange lead wire, orange button snap, orange touchproof connector

Part Number	Type	Lead Length	Quantity	Compatibility
1411517	Snap Leads	39" (1m)	2/pkg	Most systems that accept 1.5mm touchproof connectors
1411524	Snap Leads	79" (2m)	2/pkg	



### For use with Tab Style Disposable Electrodes

Part #	Item	Lead Length	Qty
301633	Alligator Clip Leads - (R/WT/BL)	18 in	3/pkg
301636	Alligator Clip Leads - (R/G/BK)	1 m	3/pkg
301635	Alligator Clip Leads - (R/WT/BL/G/BK)	18 in	5/pkg



## Disposable Adhesive Center Snap Round Electrode

- Circular shape, 1-7/8" (48mm) diameter
- Solid gel Ag/AgCl sensor material with vinyl backing


Part Number	Type	Quantity	Compatibility
022107	Circular with center snap	5/pkg	Most systems that accept 1.5mm touchproof connectors; requires snap leads
1410500	Circular with center snap	30/pkg	



## Gold TIPtrodes/Reusable TIPtrode Cable

- Gold foil-wrapped polyurethane foam eartip electrodes offer both sound delivery and recording functions in one unit
- Enhances Wave I amplitude but may be used for ABR and all other auditory EP recording/testing
- TIPtrodes require one pair of TIPtrode cables
- Single patient use only

Part Number	Type	Quantity	Compatibility
019-746802	10mm TIPtrode	20/pkg	All Natus Tubal Insert Phones
019-741901	13mm TIPtrode	20/pkg	
085-415000	9.8" (250mm) TIPtrode Cable	1 pair/pkg	



## US and Canada phone orders

**Toll Free:** 1-800-356-0007

Monday through Friday

8 am to 5 pm (CST)

## Fax or e-mail orders

24 hours a day, 7 days a week

**Fax:** 608-829-8709

**Email:** [supplies.nc@natus.com](mailto:supplies.nc@natus.com)

## US Online Orders

Natus Neuro Store

[www.natusneurostore.com](http://www.natusneurostore.com)

## International orders

**Phone:** +1 608-829-8500

**Fax:** +1 608-829-8770

**Email:** [GMB-IntlOrdersNC@natus.com](mailto:GMB-IntlOrdersNC@natus.com)

## Mail in orders

### **Natus Medical Incorporated**

P.O. Box 44994

Madison, WI 53744-4994

Attention: Customer Service Department



#### **Natus Neuro**

3150 Pleasant View Road  
Middleton, WI 53562 USA

Tel: 1-800-356-0007

1-608-829-8500

Fax: 1-608-829-8709

[www.natus.com](http://www.natus.com)

©2018 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.