



Natus EP Supplies

Insert earphones, Headphones, Bone Vibrator

AUDITORY

Tubal Insert Phones (TIPs) Kits

- ✓ For auditory studies including evoked potentials, audiometry & intraoperative monitoring
- ✓ TIPs maximize interaural stimulus isolation, help prevent ear canal collapse & reduce stimulus artifact



Reusable Left & Right Tubal Insert Phones with Cable

- Three red and three blue silicone tubes
- Two clips with Velcro loop pads
- Disposable foam ear tips:
 - Fifty regular/13mm
 - Four jumbo/18mm
 - Ten baby/10mm
- Impedance probe adapters

Part Number	Description	Quantity	Compatibility
041-704000	TIP 300 Kit	1 kit/pkg	Nicolet systems

Part Number	Description	Cable Length	Quantity
085-746400	Replacement Cable	16' (4.88m)	1/pkg



Reusable Left & Right Tubal Insert Phones with Cable

- Three red and three blue silicone tubes
- Two clips with Velcro loop pads
- Disposable foam ear tips:
 - Fifty regular/13mm
 - Four jumbo/18mm
 - Ten baby/10mm
 - Ten infant/3.5mm
 - Ten infant/4.0mm
- Impedance probe adapters

Part Number	Description	Connector	Quantity	Compatibility
10250	Xltek XCalibur EMG & Protektor IOM TIPs Kit	Lemo	1 kit/pkg	Xltek XCalibur & Protektor systems



Reusable Left & Right Tubal Insert Phones with Cable

Part Number	Description	Quantity	Compatibility
9031E0277	Dantec Keypoint TIPs Kit	1 kit/pkg	Dantec Keypoint G4, Keypoint Focus & Natus UltraPro S100 EMG systems

AUDITORY

Tubal Insert Phones (TIPs) Supplies

Disposable Polyurethane Foam Ear Tips

- Single patient use only


Part Number	Type	Quantity	Compatibility
103342	18mm, yellow	24/pkg	All Natus Tubal Insert Phones
1700-9604	13mm, yellow	50/pkg	
1700-9605	10mm, beige	50/pkg	



Replacement Red/Blue Silicone Tubes

- Three red and three blue tubes per pack
- Three pairs of additional replacement connectors
- Nubbies are designed to work with replacement tubes

Part Number	Type	Quantity	Compatibility
1700-9688	Red/Blue Tubes	3 each color/pkg	All Natus Tubal Insert Phones
1700-9686	Nubbies	10/pkg	



Replacement Clear Silicone Tubes

- Three pairs clear tubes per pack
- Three pair replacement connectors
- Nubbies are designed to work with replacement tubes

Part Number	Type	Quantity	Compatibility
019-411200	Clear Tubes	6/pkg	All Natus Tubal Insert Phones
1700-9686	Nubbies	10/pkg	



Replacement Impedance Adapter Tip


Part Number	Type	Quantity	Compatibility
019-428400	Adapter Tip	4/pkg	All Natus Tubal Insert Phones



Gold TIPtrodes/Reusable TIPtrode Cable

- Gold foil-wrapped polyurethane foam eartip electrodes offer both sound delivery and recording functions in one unit
- Enhances Wave I amplitude but may be used for ABR and all other auditory EP recording/testing
- TIPtrodes require one pair of TIPtrode cables
- Single patient use only

Part Number	Type	Quantity	Compatibility
019-746802	10mm TIPtrode	20/pkg	All Natus Tubal Insert Phones
019-741901	13mm TIPtrode	20/pkg	
085-415000	9.8" (250mm) TIPtrode Cable	1 pair/pkg	



EP HEADPHONES

Headphones

Nicolet® Headphones

Part Number	Description	Compatibility
842-202700	Unshielded Headphones, 300 Ohm	VikingQuest or EDX system with auditory stimulator
842-202300	Shielded Headphones, 300 Ohm	



842-202700



842-202300

Dantec® Keypoint® Headphones

Part Number	Description	Compatibility
9031E0258	Calibrated Headphones	Natus UltraPro S100, Dantec Keypoint Focus & Keypoint G4 with auditory stimulator





300 Ohm Bone Vibrator Transducer

- Used to determine conductive hearing loss or indicate presence of middle ear pathology
- Includes 15' (4.6m) cable

Part Number	Type	Quantity	Compatibility
842-202600	Bone Vibrator Transducer	1/pkg	Nicolet VikingQuest

US and Canada phone orders

Toll Free: 1-800-356-0007

Monday through Friday

8 am to 5 pm (CST)

Fax or e-mail orders

24 hours a day, 7 days a week

Fax: 608-829-8709

Email: supplies.nc@natus.com

US Online Orders

Natus Neuro Store

www.natusneurostore.com

International orders

Phone: +1 608-829-8500

Fax: +1 608-829-8770

Email: GMB-IntlOrdersNC@natus.com

Mail in orders

Natus Medical Incorporated

P.O. Box 44994

Madison, WI 53744-4994

Attention: Customer Service Department



Natus Neuro

3150 Pleasant View Road
Middleton, WI 53562 USA

Tel: 1-800-356-0007

1-608-829-8500

Fax: 1-608-829-8709

www.natus.com

©2018 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.